

## CO<sub>2</sub> measuring system

Both, taste and shelf life of the product are affected by the amount of CO<sub>2</sub> dissolved in the product.

A reliable measuring method is the basis for maintaining a uniform carbon dioxide content in beer, carbon-ated water and soft drinks.

The demand for an accurate and user independent CO<sub>2</sub>-measurement device led to the development of the automatic Steinfurth CO2MS.

### Description

Based on the laws of Henry and Dalton, the concentration of CO<sub>2</sub> dissolved in a liquid can be determined in a closed package by measuring the pressure and temperature when a state of equilibrium between the gas and the liquid phase exists

The CO2MS accurately computes the CO<sub>2</sub> by using these measurements in conjunction with the particular product CO<sub>2</sub>-formula.

A state of equilibrium is achieved in our process by evenly rotating the bottle in a vertical orientation for a short period of time prior to performing the measurement.

The overhead tumbling of the beverage packages provides an optimal equilibrium phase for accurate calculations of the CO<sub>2</sub> content.

The Steinfurth CO2MS combines the perfect sample preparation with accurate measurement.

### Features

- Extremely repeatable results
- User independent automatic operation
- Optimal automatic sample preparation
- Low maintenance requirements
- Easy checking and calibration
- Suitable for all packaging & beverage types
- Simulation of the beverage consuming situation
- Easy operation
- Simultaneous sample preparation and measurement
- Quick return on investment



### Technical specifications

#### Package type

Bottle or can

#### Duration of 1 measuring

approx. 1 to 2 minutes

#### Data output

RS 232 & Ethernet

#### Power supply

230 VAC / 115 VAC

#### Accuracy (pressure)

+/- 0,03 bar (0,44 PSI)

#### Accuracy (temperature)

+/-0,3 °C (0,54 °F)

#### CO<sub>2</sub> repeatability

+/- 0,05 g/l (0,025 vol)

#### Max. pressure

10 bar (145 PSI)

## Versions

### Model CO2MS- 1

The PK 1, fitted with analogue manometer, displays the equilibrium pressure at the end of the shaking process.

The temperature needs to be measured with a thermometer and the carbon dioxide content read off a CO<sub>2</sub>-chart, or calculated manually.

### Model CO2MS- 2

The PK 2, fitted with a high precision digital manometer, displays the equilibrium pressure at the end of the shaking process.

The temperature needs to be measured with a thermometer and the carbon dioxide content read off a CO<sub>2</sub>-chart, or calculated manually.

### Model CO2MS- 2C

The PK 2C, measures pressure and temperature simultaneously, calculates the CO<sub>2</sub>-content and displays all three parameters on the LCD. The last CO<sub>2</sub>-measurement result is retained in the memory even after the device is switched off and can be recalled any time.

### Model MK 6 (smart CDA)

The MK 6 measures pressure and temperature simultaneously and transfers these results to the base unit.

The measurement procedure, including preliminary shaking, head pressure release, procedure "breaks" and measurement is fully automated.

The micro controller in the base unit not only calculates and displays the carbon dioxide content on the touch screen, but also writes all the data in form of a measurement test report into the internal memory.

The measurement results content CO<sub>2</sub>-result, pressure, temperature, date, time, identification code of the measuring head and sample ID code.

One input possibility for individual customized sample ID codes is possible via barcode scanner. The measuring data can be transferred from the system via RS 232 or network interface directly into the evaluation software or customized IT system.

## Accessories

### Bottle-/ Can Insert / PET-Neck Holder



### Software / PC-Cable



### Barcode Scanner

