



rycobelgroup

Proven quality

RYC-D-09-36-M222QW

QUICKWASH PLUS



The QuickWash Plus gives fabric shrinkage information in a 12 15 minute test versus the 8 12 hours it takes a standard washer and dryer. This instrument is ideal for the factory floor and provides rapid, on the spot information on fabric in production so that the stenter frame, compactor, and sanforizer can be set to maximize fabric yield. This allows user to bypass the standard wash and dry cycles that take so long thus significantly reducing labor in the lab.

How the QuickWash saves users money

- Eliminates unnecessary production cost of overconditioning fabric
- Maximizes fabric yield
- Eliminates most standard washer and dryer cycles which reduces:
 - Labor
 - Water usage
 - Electricity
 - Detergent Usage

A chambered basket is agitated at high speed during the wash, rinse, and extraction cycles. Samples are then subjected to a controlled flow of heated air which causes a tumbling action while an infrared sensor measures the samples' surface temperature.

Technical specifications

- Width: 685 mm
- Depth: 635 mm
- Height 585 mm
- Weight: 137 kg